

# YEASTCYTE<sup>®</sup> Bottle and Tube Kit

YeastCyte<sup>®</sup> Bottle and Tube kit contain the following and are designed to be used with YeastCyte Reagent Kit.

Prod. No. 3020	Description	Used for
Container A	10000ml, plastic PE, with tube and filter	Sheat Liquid
Container B	10000ml, plastic PE, with tube	Low waste
Bottle A	1000ml, clear glass with tube	High concentrate waste
Bottle B	1000ml, brown Duran with tube and filter	Dye solution
Column*	With Dowex to prevent propidium into the waste	High concentrate waste
Dowex*	Filter material	The column

\* Optional, according to Internal requirements.

## Key advantages:

- Easy to use, all tubes , filters and corks ready with connectors to be inserted into instrument and bottles
- The set can be used with refill from YeastCyte reagent kit
- Easy to install on the table or in a suitable shelf under the table



*Bottle and Tube kit*

Every kit is carefully designed to fit with the YeastCyte Reagent kit, to ensure a safe and trouble free flow of liquids to and from the YeastCyte instrument. By using the optional column and Dowex on the high concentrated waste you will also easy take care of the environment in a safety way.

More information on [www.biodetect.no](http://www.biodetect.no) .